


A2-DX


NIPPON SIGMAX CO.,LTD. 33rd Floor Shinjuku Oak Tower 6-8-1 Nishi-shinjuku, Shinjuku-ku, Mone: +81 3 5326 3229











ZAMST A2-DX ANKLE SUPPORT INSTRUCTIONS FOR USE

Please read through the following instructions completely before using the product.

Intended use: Ankle protection. Restriction of ankle inversion and eversion. Anterior ankle stability.

- *This product is not for medical treatment.
- *This product is for ankle support; however, it will not prevent all ankle injuries.

In use of the product:

- The product may become discolored due to the storage and/or usage conditions. However, it does not affect on
 product functions.
- Do not put excess weight on the hook-and-loop fastener and the raised fabric when store the product. The hook-and-loop fastener may not work effectively.
- Do not use the product while driving vehicle. It may interfere with ability to drive.

Do not use this product in the following cases:

 If your skin is sensitive or you have allergies, you may develop a rash or allergic reactions while using the product. If such symptoms occur, discontinue use and consult with a physician.

 - If you have abnormal conditions such as cuts, fractures, humbness, swelling or rash in the area where product is used.

Safety precautions:

- Do not use the product other than ankle joint.
- Make sure to use the correct size and apply the support to the correct foot (left or right). (Refer to the size chart on the package and size tag inside the product)
- -Wear the product over a sock. Wearing the product in direct contact with skin may cause skin discomfort.
- Do not use the product while sleeping. If the product shifts to an abnormal position or is used for extended periods
 while sleeping, it may interfere with blood circulation.
- -If the body bar is compressed for extended periods, blood circulation may be adversely affected or you may experience discomfort. To avoid such conditions, it is recommended that you remove and re-apply the product at regular intervals.
- -Make sure to use the correct position; top or bottom, inside or outside, front or back.
- -Do not pull the strap too tight or create excessive compression with the support. Blood circulation may be adversely affected or you may experience skin problem. -If the product shifts out of place or becomes too loose, re-apply the product starting with the first step described in
- If the product shifts out of place or becomes too loose, re-apply the product starting with the first step described in the "How to apply" instructions.
- -If the product is broken or worn, or if the hook-and-loop fastener does not work, do not use the product.
- -Do not modify or repair the product. The product may lose its effectiveness and/or cause injuries.
 -If you play sports where the product's plastic components (such as rigid guards) could cause harm to other people, do not use the product.
 - -If you experience pain, swelling, skin irritation, discomfort or any other abnormal reaction while using this product, discontinue use and consult with a physician.
 - -The product may become discolored by perspiration or other liquids. The product' s dye may bleed onto clothing worn under or near the product.

How to wash:

- -The product's dye may bleed. Do not wash with whites or bright colors.
- -Make sure that the hook-and-loop fasteners do not touch fabric. They may damage fabric.
- -Hand wash with mild detergent and air dry in the shade.
- -Do not place in dryer. Do not dry clean or iron. Do not use products containing chlorine, including bleach. These may damage the product.

How to apply:

- *Illustrations are for the right ankle.
- Release all hook-and loop fasteners. Lay the rigid guards flat behind the ankle. Hold the ankle in a vertical position. Slide foot into the support from the back <u>until the inner ankle hone is aligned with the area indicated by the</u> <u>dashed circle</u> which is right above the rivet on the rigid guard. (Illustration 1a)
- Make sure that the center line of the foot is positioned just under the support's center line (seam). (Illustrations 1b)
- 2. Pull hook A and hook tab B upward at opposite angles. Fasten hook A. (Illustration 2)
- 3. Wrap the support tightly around the ankle by pulling hook B and hook fastener C so as not to loosen the support. Fasten B. (Illustration 3)
- Pull the crossing straps at opposing angles. Fasten their edges above the ankle bone on opposite sides of the ankle. (Illustration 4)
- 5. Stand the rigid guards, upright. Position them so that the raised area on the guard (just above the rivet) comes over the bone of the ankle which is indicated by the dashed circle D, by adjusting the length of the rear position of the upper strap. Fasten the hooks located on the inner side of each rigid guard to stabilize the guards. (Illustration 5a & 5b)
- 6. Use one hand to prevent the rigid guards from moving. Wrap the upper strap around the leg and pass the strap through the plastic loop E. (Illustration 6a) Fold the strap back over itself and fasten hook. Do not pull the strap too tight or create excessive compression. (Illustration 6b)
- Product properly fastened. Illustration a. shows logo in the center, however, the logo position does not guide the proper fastening position of the upper strap. It does not always come in the center and different with each individual. (Illustration 7)
- a. Front view b. Lateral view c. Back view
 8. Inserting buffers F.
 - This product comes with two buffer accessories F. If you experience discomfort caused by the rigid materials, use the buffers. (Illustration 8)
- Example: a. For top part of the rigid guard
- b. For ankle bone area
- c. For inside rivet area

Composition: Fiber: Nylon, Polyester, Elastane

Non-Fiber: Chloroprene, Nylon, Polyurethane, Polyethylene, Polyacetal, Imitation leather, Brass