




TS-1


NIPPON SIGMAX CO., LTD.
 33rd Floor Shinjuku Oak Tower
 6-8-1 Nishi-shinjuku, Shinjuku-ku,
 Tokyo 163-6033 JAPAN
 Phone: +81 3 5326 3229




2-a




2-b




3



4



5



6

ZAMST TS-1 THIGH SUPPORT INSTRUCTIONS FOR USE

Please read through the following instructions completely before using the product.

Intended use: Thigh compression

*This product is not for medical treatment.

*This product is for thigh support; however, it will not prevent all thigh injuries.

Do not use this product in the following cases:

- If you have abnormal conditions such as cuts, fractures, numbness, swelling or rash in the area where product is used.

Safety precautions:

- If you experience pain, swelling, skin irritation, discomfort or any other abnormal reaction while using this product, discontinue use and consult with a physician.

- If your skin is sensitive or you have allergies, you may develop a rash or allergic reactions while using the product. If such symptoms occur, discontinue use and consult with a physician.

- Do not use the product while sleeping. If the product shifts to an abnormal position or is used for extended periods while sleeping, it may interfere with blood circulation.

- If the body part is compressed for extended periods, blood circulation may be adversely affected or you may experience discomfort. To avoid such conditions, it is recommended that you remove and re-apply the product at regular intervals.

- If the product is broken or worn or if the hook-and-loop fastener does not work, do not use the product.

- Make sure that you use the correct size. (Refer to the size chart on the package and size tag inside the product.)

- Do not pull the straps too tight or create excessive compression with the support.

- If the product shifts out of place or becomes too loose, re-apply the product starting with the first step described in the "How to apply" instructions.

- Do not modify or repair the product. The product may lose its effectiveness and/or cause injuries.

To prevent product damage:

- The product may become discolored by perspiration or other liquids. The product's dye may bleed onto clothing worn under or near the product.

- Make sure that the hook-and-loop fasteners do not touch fabric. They may damage fabric.

How to apply:

* The product is not designed specifically for use on the right or left, so application is the same for both legs.

1. Release all hook and loop fasteners. (Illustration 1)

2. Place the main body (A in Illustration 1) on the area where you are experiencing discomfort. Keep the seam parallel to the leg.

Illustration 2-a: Front thigh

Illustration 2-b: Back thigh

3. Pull both straps around the leg. Make sure that the main body does not twist or shift away from the area of discomfort. (In Illustration 3 the area of discomfort is at the back of the thigh.)

4. Pull both straps horizontally in opposite directions. (Illustration 4)

5. As you continue pulling, fasten the hook on the lower strap. (Illustration 5)

6. As you continue pulling, fasten the hook on the upper strap. (Illustration 6)

How to wash:

- Hand wash with mild detergent and air dry.

- Do not place in dryer. Do not dry clean or iron. Do not use products containing chlorine, including bleach. These may damage the product.

- The product's dye may bleed. Do not wash with whites or bright colors.

Composition: Polystyrene, Nylon, Polyester, Polyurethane